



Exam Factory Documentation

Overview

The **Exam Factory** (`exam_factory.py`) is the core utility that **builds an Exam object from a JSON template**.

It ties together **backend models** (`Exam`, `Segment`, `Measurement`) with **template definitions** (`carotid.json`, `renal.json`, etc.).

This function is used whenever a **new ultrasound exam is created** in the Lumen system — whether it comes from: - **EPIC HL7 arrival messages** (patient scheduled and marked “arrived”)

- **Manual entry by a technologist**

- **Agent-based test cases or imports**



Responsibilities

The factory has **one clear responsibility**:

Given an exam type and site, generate a new `Exam` record with all its segments and measurement placeholders initialized according to the template.



Workflow (Step-by-Step)

1. Load Template

```
template = get_template(exam_type, site)
```

- Uses `template_registry` to load the correct JSON (e.g., `carotid.json`).

- Template defines: - Segments (prox ICA, mid ICA, vertebral, etc.) - Which measurements each segment has (PSV, EDV, ICA/CCA ratio, etc.) - Units for each measurement - Display groups (Right Side, Left Side, Temporal Arteries)

2. Create Exam

```
exam = Exam.objects.create(
    patient_name=patient_name,
    gender=gender,
    mrn=patient_data.get("mrn", ""),
```

```

        exam_type=exam_type,
        ...
    )

```

- Creates the **parent Exam record** in the database.
- Stores patient metadata, scope, extent, CPT, ICD-10 code, etc.
- Sets initial status = "draft".
- Patient demographics can come from **EPIC HL7 feeds** or manual entry.

3. Create Segments

```

segment = Segment.objects.create(
    exam=exam,
    name=seg["id"],
    artery=seg["vessel"].lower(),
    side=seg.get("side", "n/a")
)

```

- Each **anatomical segment** (prox ICA, mid ICA, etc.) becomes a `Segment` DB row.
- Linked to the parent `Exam`.
- `side` is stored as "left", "right", or "temporal".

4. Create Measurements

```

measurement = Measurement.objects.create(segment=segment)

```

- For each segment, a **Measurement object** is initialized.
- Defines all the numerical values, categorical flags, and dropdowns.
- Examples: - PSV (cm/s)
- EDV (cm/s)
- ICA/CCA ratio (ratio)
- Plaque morphology
- Direction of flow

5. Initialize Fields

```

for m in seg.get("measurements", []):
    field_name = m if isinstance(m, str) else m.get("name")

    if hasattr(measurement, field_name):

```

```

# Core field (DB column)
setattr(measurement, field_name, None)
else:
    # Non-core field → store in JSON
    measurement.additional_data[field_name] = None

```

- Supports two styles of template definition:
- **Compact:** `["psv", "edv"]`
- **Verbose:** `[{"name": "psv", "unit": "cm/s"}]`
- Core fields (psv, edv, ica_cca_ratio) → saved directly as DB columns.
- Non-core (artery_diameter, ap_tr, longitudinal, etc.) → saved in `additional_data` JSON field.

6. Store Units

```

if field_unit:
    measurement.additional_data[f"{field_name}_unit"] = field_unit

```

- Every measurement has a **unit** (cm/s, ratio, cm, etc.).
- Units stored in `additional_data`, so they are always available for rendering in UI or PDF.

Data Model Relationships

Here's the **mental model** for what the factory builds:

```

Exam (Carotid Exam for John Doe)
├── Segment (Right ICA Proximal)
│   └── Measurement (PSV=None, EDV=None, Ratio=None, additional_data={})
├── Segment (Right ICA Mid)
│   └── Measurement (PSV=None, EDV=None, Ratio=None)
├── Segment (Left ICA Proximal)
│   └── Measurement (PSV=None, EDV=None, Ratio=None)
└── ...

```

How It's Used

- **When HL7 says "Exam Arrived":**

```
create_exam_from_template("carotid", "mount_sinai_hospital", patient_data, tech_user)
```

- **When tech manually creates exam:**

Same call, but `patient_data` comes from frontend form.

- **When testing with sample patients:**

Call factory directly in scripts or management commands.

Example Call

```
exam = create_exam_from_template(  
    exam_type="carotid",  
    site="mount_sinai_hospital",  
    patient_data={  
        "name": "Jane Doe",  
        "gender": "female",  
        "mrn": "123456",  
        "dob": "1970-01-01",  
        "accession": "ACC-98765",  
        "scope": "bilateral",  
        "extent": "complete",  
        "cpt_code": "93880",  
        "technique": "Duplex carotid ultrasound",  
        "operative_history": "Prior left CEA",  
        "indication": "I65.23" # ICD-10 code  
    },  
    created_by="tech_001"  
)
```

Result in DB: - Exam record created

- ~20 Segment rows (prox/mid/dist ICA, CCA, vertebral, subclavian, etc.)

- Each with 1 Measurement row initialized with empty values

Key Takeaways for New Engineers

- The **factory is the single point of truth** for creating a new exam.
- It ensures template-driven **consistency**: same segments, same fields every time.
- **Core fields** = DB columns.

- **Non-core fields** = `additional_data` JSON.
- **Units** always stored alongside values.
- This enables:
 - Flexible frontends (auto-build forms from JSON)
 - Reliable backends (calculators and reports always have the right fields)
 - Easy onboarding for new exam types (just add a template JSON)